

SPECIFICATIONS
DOC 310-PR2336797

Page 1 of 2

ITEM NO.	DESCRIPTION	CHECK ✓ IF AS SPECIFIED	DESCRIBE FULLY IF NOT AS SPECIFIED
1.	Gurney, Stryker Power-Pro XT, P/N: 6506 or equal		
1.1	Brand New		
1.2	Power load compatible		
1.4	LENGTH: 82" Maximum – 62" Minimum		
1.5	WIDTH: No less than 23"		
1.6	HEIGHT RANGE: Minimum 14" to a Maximum of 42". Ability to adjust and position at any height between minimum and maximum.		
1.7	Backrest shall change positions from 0 to 73 degrees.		
1.8	WIDTH: 24" MAX - 22" MIN		
1.9	WEIGHT: Maximum gurney weight of 135 lbs. without mattress		
1.10	All Power-Pro XT standard components included		
1.11	Gurney shall be fully operational and safe while carrying a maximum weight of 700 lbs.		
1.12	Wheels shall not be less than 6" diameter by 2" wide.		
1.13	Dual wheel lock		
1.14	Standard fowler		
1.15	Patient right cot retaining post		
1.16	No runner		
1.17	Retractable head section oxygen bottle holder		
1.18	Head end storage flat		
1.19	J Hook		
1.20	X-Restraint package		
1.21	Patient right 3-Stage IV Pole		
1.22	Equipment hook		

SPECIFICATIONS DOC 310-PR2336797

Page 2 of 2

1.23	XPS Side Rails		
1.24	XPS Bolster Mattress		
1.25	Trendelenburg Option		
1.26	Dual Compatibility		
1.27	Backrest storage pouch		
1.28	120V AC, 12V DC smart charging kit		
1.29	UL Listed		
1.30	User manual, written in English.		
1.31	Warranties include travel, parts & labor		
1.32	2-Year bumper to bumper warranty		
1.33	3-year X-Frame Powertrain Warranty		
1.34	Assembly/set-up included		
2.	Shipping/Delivery		
2.1	Deliver Two (2) cots M-F between 7:30 & 2:30 pm to:		
	Monroe Correctional Complex 16500 177 th Ave SE Monroe, WA 99362		
2.2	Deliver One (3) cot M-F between 7:30 & 3:00 pm to:		
	Stafford Creek Correction Center 191 Constantine Way Aberdeen, WA 98520		
2.3	Deliver One (2) cot M-F between 7:30 & 3:00 pm to:		
	Washington State Penitentiary Attn: Warehouse 1313 N. 13 th Avenue Walla Walla, WA 99362		
3.	Delivery of Goods		
3.1	Gurneys must be delivered to each facility no later than June 30, 2024.		